


Vector Summary

Vector ID	VB190928-1144jva			
Vector Name	pPB[Exp]-Bsd-EFS>hEEF2K[NM_013302.4]*			
Date Created (Pacific Time)	2019-09-28			
Vector Size	8088 bp			
Vector Type	Mammalian Gene Expression PiggyBac Vector			
Inserted Promoter	EFS			
Inserted ORF	hEEF2K[NM_013302.4]*			
Inserted Marker	Bsd			
Plasmid Copy Number	High			
Antibiotic Resistance	Ampicillin			
Cloning Host	VB UltraStable (or alternative strain)			

Vector Map



Vector Components

Vector #: VB190928-1144jva Page 1/5



Name	Position	Size (bp)	Туре	Description	Application notes
5' ITR	■ 1-313	313	ITR	piggyBac 5' inverted terminal repeat	Recognized by PBase transposase; DNA flanked by piggyBac 5' ITR and 3' ITR can be transposed by PBase into TTAA sites.
EFS	■ 337-568	232	Promoter	Human eukaryotic translation elongation factor 1 α1 short form	Strong promoter.
Kozak	593-598	6	Miscellaneous	Kozak translation initiation sequence	Facilitates translation initiation of ATG start codon downstream of the Kozak sequence.
hEEF2K[NM_013302.4]*	■ 599-2776	2178	ORF	None	None
rBG pA	■ 2804-3325	522	PolyA_signal	Rabbit beta-globin polyadenylation signal	Allows transcription termination and polyadenylation of mRNA transcribed by Pol II RNA polymerase.
CMV promoter	■ 3351-3938	588	Promoter	Human cytomegalovirus immediate early enhancer/promoter	Strong promoter; may have variable strength in some cell types.
Bsd	■ 3970-4368	399	ORF	Blasticidin resistance gene	Allows cells to be resistant to blasticidin.
BGH pA	4 412-4636	225	PolyA_signal	Bovine growth hormone polyadenylation signal	Allows transcription termination and polyadenylation of mRNA transcribed by Pol II RNA polymerase.
3' ITR	complement (4818- 5052)	235	ITR	piggyBac 3' inverted terminal repeat	Recognized by PBase transposase; DNA flanked by piggyBac 5' ITR and 3' ITR can be transposed by PBase into TTAA sites.
Ampicillin	5884-6744	861	ORF	Ampicillin resistance gene	Allows E. coli to be resistant to ampicillin.

Vector #: VB190928-1144jva Page 2 / 5



Name	Position	Size (bp)	Туре	Description	Application notes
pUC ori	■ 6915-7503	589	Rep_origin	pUC origin of replication	Facilitates plasmid replication in E. coli; regulates high-copy plasmid number (500- 700).

Note: Components added by user are listed in **bold red** text.

Vector Sequence

1	TTAACCCTAG	AAAGATAGTC	TGCGTAAAAT	TGACGCATGC	ATTCTTGAAA	TATTGCTCTC	TCTTTCTAAA	TAGCGCGAZ
101	CATCTCAGTC	GCCGCTTGGA	GCTCCCGTGA	GGCGTGCTTG	TCAATGCGGT	AAGTGTCACT	GATTTTGAAC	TATAACGA(
201	TGATTATCTT	TTACGTGACT	TTTAAGATTT	AACTCATACG	ATAATTATAT	TGTTATTTCA	TGTTCTACTT	ACGTGATA
301	TGTTATAGAT	ATCATCAACT	TTGTATAGAA	AAGTTGGGCT	CCGGTGCCCG	TCAGTGGGCA	GAGCGCACAT	CGCCCACA
401	GGGTCGGCAA	TTGATCCGGT	GCCTAGAGAA	GGTGGCGCGG	GGTAAACTGG	GAAAGTGATG	TCGTGTACTG	GCTCCGCC'
501	ACCGTATATA	AGTGCAGTAG	TCGCCGTGAA	CGTTCTTTTT	CGCAACGGGT	TTGCCGCCAG	AACACAGGCA	AGTTTGTA(
601	GGCAGACGAA	GATCTCATCT	TCCGCCTGGA	AGGCGTTGAT	GGCGGCCAGT	CCCCCGAGC	TGGCCATGAT	GGTGATTC!
701	GAAGGTTACT	TCATCTGCCC	CATCACGGAT	GACCCAAGCT	CGAACCAGAA	TGTCAATTCC	AAGGTTAATA	AGTACTAC?
801	GGTATAGCTC	CAGCGGGTCC	CCGGCAAACT	CCTTCCACTT	CAAGGAAGCC	TGGAAGCACG	CAATCCAGAA	GGCCAAGC?
901	GTTCCACCTG	GAAGATATTG	CCACCGAACG	TGCTACTCGA	CACAGGTACA	ACGCCGTCAC	CGGGGAATGG	CTGGATGA!
1001	TCTCAGCCCT	TCGGCCGAGG	AGCAATGAGG	GAGTGCTTCC	GGACGAAGAA	GCTCTCCAAC	TTCTTGCATG	CCCAGCAG!
1101	TGGCGAAGCG	CTACATCGAG	CCCGTAGACC	GGGATGTGTA	CTTTGAGGAC	GTGCGTCTAC	AGATGGAGGC	CAAGCTCT(
1201	CAAGCCCCCC	AAGCAGGTGG	ACATCATGCA	GATGTGCATC	ATCGAGCTGA	AGGACAGACC	GGGCAAGCCC	CTCTTCCA
1301	AAGTACATCA	AGTACAACTC	CAACTCTGGC	TTTGTCCGCG	ATGACAACAT	CCGCCTGACG	CCGCAGGCCT	TCAGCCAC'
1401	ATCAGCTGAT	AGTGGTGGAC	ATCCAGGGAG	TTGGGGATCT	CTACACTGAC	CCACAGATCC	ACACGGAGAC	GGGCACTG2
1501	TGTCCGCGGG	ATGGCGCTCT	TCTTCTACTC	TCATGCCTGC	AACCGGATTT	GCGAGAGCAT	GGGCCTTGCT	CCCTTTGA(
1601	GCAGTGAATC	AGAACACCAA	GCTGCTGCAA	TCAGCCAAGA	CCATCTTGAG	AGGAACAGAG	GAAAAATGTG	GGAGCCCC
1701	GCCGGCCACC	CCTGCTCCGT	CCCCTTTCAG	AGAACTCTGG	AGACGAGAAC	ATGAGCGACG	TGACCTTCGA	CTCTCTCC
1801	ACCACACAGC	CAGAAGCTAG	ACCACCTCCA	TTGGCCAGTG	TTCAGTGACC	TCGATAACAT	GGCATCCAGA	GACCATGA'
1901	TCTGAGAATA	GTGGGGACAG	CGGATACCCC	AGTGAGAAGC	GGGGTGAGCT	GGATGACCCT	GAGCCCCGAG	AACATGGC
2001	ACGAGTCTGA	CGAAGACAGC	CTGGGCAGCT	CTGGACGGGT	ATGTGTAGAG	AAGTGGAATC	TCCTCAACTC	CTCCCGCC'.
2101	CGTGGCCCTG	GAAGTGCAAA	GGCTTAATGC	TCTGGACCTC	GAAAAGAAAA	TCGGGAAGTC	CATTTTGGGG	AAGGTCCA'
2201	GAGGGTGGGC	GCTTCTGCGA	GAAGGGCGAG	GAGTGGGACC	AGGAGTCGGC	TGTCTTCCAC	CTGGAGCACG	CAGCCAAC
2301	TGGGCCTGGG	ACTCATGTAC	TCGCAGTTGC	CTCATCACAT	CCTAGCCGAT	GTCTCTCTGA	AGGAGACAGA	AGAGAACA
2401	ACTAAAGGCC	GCTGAAGCTG	GCGACAGGCA	GTCCATGATC	CTAGTGGCGC	GAGCTTTTGA	CTCTGGCCAG	AACCTCAG
2501	CTAGAGGCCC	TGCACTGGTA	CAACACTGCC	CTGGAGATGA	CGGACTGTGA	TGAGGGCGGT	GAGTACGACG	GAATGCAG
2601	TGGCCAGGGA	GGCCGAGATG	CTGTTCACAG	GAGGCTACGG	GCTGGAGAAG	GACCCGCAGA	GATCAGGGGA	CTTGTATA(
2701	GGAAGCCATG	AAGGGCCGAC	TGGCCAACCA	GTACTACCAA	AAGGCTGAAG	AGGCCTGGGC	CCAGATGGAG	GAGTAAAC(
2801	TGATCCTCAG	GTGCAGGCTG	CCTATCAGAA	GGTGGTGGCT	GGTGTGGCCA	ATGCCCTGGC	TCACAAATAC	CACTGAGA:
2901	TATGGGGACA	TCATGAAGCC	CCTTGAGCAT	CTGACTTCTG	GCTAATAAAG	GAAATTTATT	TTCATTGCAA	TAGTGTGT'
3001	TCGGAAGGAC	ATATGGGAGG	GCAAATCATT	TAAAACATCA	GAATGAGTAT	TTGGTTTAGA	GTTTGGCAAC	ATATGCCC2

Vector #: VB190928-1144jva Page 3 / 5



3101	AGGTTGGCTA	TAAAGAGGTC	ATCAGTATAT	GAAACAGCCC	CCTGCTGTCC	ATTCCTTATT	CCATAGAAAA	GCCTTGAC:
3201	TTTTGTTTTG	TGTTATTTTT	TTCTTTAACA	TCCCTAAAAT	TTTCCTTACA	TGTTTTACTA	GCCAGATTTT	TCCTCCTC'
3301	CTGTCCCTCT	TCTCTTATGG	AGATCCCTCG	ACCTGCAGCC	CAAGCTTCGC	GTTGACATTG	ATTATTGACT	AGTTATTA
3401	TTAGTTCATA	GCCCATATAT	GGAGTTCCGC	GTTACATAAC	TTACGGTAAA	TGGCCCGCCT	GGCTGACCGC	CCAACGAC(
3501	TGACGTATGT	TCCCATAGTA	ACGCCAATAG	GGACTTTCCA	TTGACGTCAA	TGGGTGGAGT	ATTTACGGTA	AACTGCCC2
3601	TCATATGCCA	AGTACGCCCC	CTATTGACGT	CAATGACGGT	AAATGGCCCG	CCTGGCATTA	TGCCCAGTAC	ATGACCTTA
3701	TACATCTACG	TATTAGTCAT	CGCTATTACC	ATGGTGATGC	GGTTTTGGCA	GTACATCAAT	GGGCGTGGAT	AGCGGTTT(
3801	TCCACCCCAT	TGACGTCAAT	GGGAGTTTGT	TTTGGCACCA	AAATCAACGG	GACTTTCCAA	AATGTCGTAA	CAACTCCG
3901	TAGGCGTGTA	CGGTGGGAGG	TCTATATAAG	CAGAGCTCTC	TGGCTAACTA	GAGAACCCAC	TGCGCCACCA	TGGCCAAG(
4001	CCCTCATTGA	AAGAGCAACG	GCTACAATCA	ACAGCATCCC	CATCTCTGAA	GACTACAGCG	TCGCCAGCGC	AGCTCTCT(
4101	TGGTGTCAAT	GTATATCATT	TTACTGGGGG	ACCTTGTGCA	GAACTCGTGG	TGCTGGGCAC	TGCTGCTGCT	GCGGCAGC'
4201	GCGATCGGAA	ATGAGAACAG	GGGCATCTTG	AGCCCCTGCG	GACGGTGCCG	ACAGGTGCTT	$\underline{\mathtt{CTCGATCTGC}}$	ATCCTGGG2
4301	GTGATGGACA	GCCGACGGCA	$\underline{GTTGGGATTC}$	GTGAATTGCT	$\underline{\mathtt{GCCCTCTGGT}}$	TATGTGTGGG	AGGGCTAACT	CGAGTCTA(
4401	ATCAGCCTCG	ACTGTGCCTT	CTAGTTGCCA	GCCATCTGTT	$\underline{\mathtt{GTTTGCCCCT}}$	$\underline{\mathtt{CCCCGTGCC}}$	$\underline{\mathtt{TTCCTTGACC}}$	CTGGAAGG'
4501	TAATAAAATG	AGGAAATTGC	ATCGCATTGT	CTGAGTAGGT	$\underline{\mathtt{GTCATTCTAT}}$	TCTGGGGGGT	GGGGTGGGGC	AGGACAGC2
4601	ATAGCAGGCA	TGCTGGGGAT	GCGGTGGGCT	CTATGGCTCG	AGTTAATTAA	CGAGAGCATA	ATATTGATAT	GTGCCAAA(
4701	GTATAATTTG	TTTCTATTAT	GTATAGGTTA	AGCTAATTAC	TTATTTTATA	ATACAACATG	ACTGTTTTTA	AAGTACAA
4801	AGAGAATGTT	TAAAAGTTTT	GTTACTTTAT	AGAAGAAATT	TTGAGTTTTT	GTTTTTTTT	AATAAATAAA	TAAACATA
4901	TTATTAGTAT	GTAAGTGTAA	ATATAATAAA	ACTTAATATC	TATTCAAATT	AATAAATAAA	CCTCGATATA	CAGACCGA!
5001	CGCATGATTA	TCTTTAACGT	ACGTCACAAT	ATGATTATCT	TTCTAGGGTT	AAATAATAGT	TTCTAATTTT	TTTATTATT
5101	GAGCTCCAAT	TCGCCCTATA	GTGAGTCGTA	TTACAATTCA	CTGGCCGTCG	TTTTACAACG	TCGTGACTGG	GAAAACCC'
5201	CTTGCAGCAC	ATCCCCCTTT	CGCCAGCTGG	CGTAATAGCG	AAGAGGCCCG	CACCGATCGC	CCTTCCCAAC	AGTTGCGC2
5301	CGCCCTGTAG	CGGCGCATTA	AGCGCGGCGG	GTGTGGTGGT	TACGCGCAGC	GTGACCGCTA	CACTTGCCAG	CGCCCTAG(
5401	CCCTTCCTTT	CTCGCCACGT	TCGCCGGCTT	TCCCCGTCAA	GCTCTAAATC	GGGGGCTCCC	TTTAGGGTTC	CGATTTAG!
5501	AAAAAACTTG	ATTAGGGTGA	TGGTTCACGT	AGTGGGCCAT	CGCCCTGATA	GACGGTTTTT	CGCCCTTTGA	CGTTGGAG!
5601	TCTTGTTCCA	AACTGGAACA	ACACTCAACC	CTATCTCGGT	CTATTCTTTT	GATTTATAAG	GGATTTTGCC	GATTTCGG(
5701	GATTTAACAA	AAATTTAACG	CGAATTTTAA	CAAAATATTA	ACGCTTACAA	TTTAGGTGGC	ACTTTTCGGG	GAAATGTG(
5801	TTTTCTAAAT	ACATTCAAAT	ATGTATCCGC	TCATGAGACA	ATAACCCTGA	TAAATGCTTC	AATAATATTG	AAAAAGGAA
5901	CCGTGTCGCC	CTTATTCCCT	TTTTTGCGGC	ATTTTGCCTT	CCTGTTTTTG	CTCACCCAGA	AACGCTGGTG	AAAGTAAA
6001	GCACGAGTGG	GTTACATCGA	ACTGGATCTC	AACAGCGGTA	AGATCCTTGA	GAGTTTTCGC	CCCGAAGAAC	GTTTTCCA
6101	TGCTATGTGG	CGCGGTATTA	TCCCGTATTG	ACGCCGGGCA	AGAGCAACTC	GGTCGCCGCA	TACACTATTC	TCAGAATG
6201	CACAGAAAAG	CATCTTACGG	ATGGCATGAC	AGTAAGAGAA	TTATGCAGTG	CTGCCATAAC	CATGAGTGAT	AACACTGC
6301	ATCGGAGGAC	CGAAGGAGCT	AACCGCTTTT	TTGCACAACA	TGGGGGATCA	TGTAACTCGC	CTTGATCGTT	GGGAACCG
6401	ACGACGAGCG	TGACACCACG	ATGCCTGTAG	CAATGGCAAC	AACGTTGCGC	AAACTATTAA	CTGGCGAACT	ACTTACTC!
6501	AGACTGGATG	GAGGCGGATA	AAGTTGCAGG	ACCACTTCTG	CGCTCGGCCC	TTCCGGCTGG	CTGGTTTATT	GCTGATAA
6601							TACACGACGG	
6701	GACAGATCGC	TGAGATAGGT	GCCTCACTGA	TTAAGCATTG	GTAACTGTCA	GACCAAGTTT	ACTCATATAT	ACTTTAGA!
6801	ATTTAAAAGG	ATCTAGGTGA	AGATCCTTTT	TGATAATCTC	ATGACCAAAA	TCCCTTAACG	TGAGTTTTCG	TTCCACTG
6901	ATCAAAGGAT	CTTCTTGAGA	$\frac{\texttt{TCCTTTTTTT}}{}$	CTGCGCGTAA	TCTGCTGCTT	GCAAACAAAA	AAACCACCGC	TACCAGCG
7001	AGCTACCAAC	TCTTTTTCCG	AAGGTAACTG	GCTTCAGCAG	AGCGCAGATA	CCAAATACTG	TTCTTCTAGT	GTAGCCGT
7101	CTCTGTAGCA	CCGCCTACAT	ACCTCGCTCT	GCTAATCCTG	TTACCAGTGG	CTGCTGCCAG	TGGCGATAAG	TCGTGTCT'

Vector #: VB190928-1144jva Page 4 / 5



7201	TAGTTACCGG	ATAAGGCGCA	GCGGTCGGGC	TGAACGGGGG	GTTCGTGCAC	ACAGCCCAGC	TTGGAGCGAA	CGACCTAC
7301	GTGAGCTATG	AGAAAGCGCC	ACGCTTCCCG	AAGGGAGAAA	GGCGGACAGG	TATCCGGTAA	GCGGCAGGGT	CGGAACAG
7401	AGGGGGAAAC	GCCTGGTATC	TTTATAGTCC	TGTCGGGTTT	CGCCACCTCT	GACTTGAGCG	TCGATTTTTG	TGATGCTC
7501	AAAAACGCCA	GCAACGCGGC	CTTTTTACGG	TTCCTGGCCT	TTTGCTGGCC	TTTTGCTCAC	ATGTTCTTTC	CTGCGTTA!
7601	GTATTACCGC	CTTTGAGTGA	GCTGATACCG	CTCGCCGCAG	CCGAACGACC	GAGCGCAGCG	AGTCAGTGAG	CGAGGAAG(
7701	ACCGCCTCTC	CCCGCGCGTT	GGCCGATTCA	TTAATGCAGC	TGGCACGACA	GGTTTCCCGA	CTGGAAAGCG	GGCAGTGA
7801	TAGCTCACTC	ATTAGGCACC	CCAGGCTTTA	CACTTTATGC	TTCCGGCTCG	TATGTTGTGT	GGAATTGTGA	GCGGATAA(
7901	ATGACCATGA	TTACGCCAAG	CTCGAAATTA	ACCCTCACTA	AAGGGAACAA	AAGCTGGTAC	CTCGCGCGAC	TTGGTTTG(
8001	TCACACAGCT	TGGCCACAAT	GTGGTTTTTG	TCAAACGAAG	ATTCTATGAC	GTGTTTAAAG	TTTAGGTCGA	GTAAAGCG(

Deleted from original sequence:

1087-3265

C

Validation by Restriction Enzyme Digestion

Restriction Enzymes	Cutting Sites	DNA Fragments (bp)
NaeI	1704, 5426	3722, 4366
NdeI	3012, 3072, 3080, 3604	60, 8, 524, 7496
DrdI	106, 4058, 5576, 7458	3952, 1518, 1882, 736
ApaLI	6000, 7246	1246, 6842
ApaLI+NaeI	1704, 5426, 6000, 7246	3722, 574, 1246, 2546
ApaLI+NdeI	3012, 3072, 3080, 3604, 6000, 7246	60, 8, 524, 2396, 1246, 3854
ApaLI+DrdI	106, 4058, 5576, 6000, 7246, 7458	3952, 1518, 424, 1246, 212, 736

Vector #: VB190928-1144jva Page 5 / 5