


Vector Summary

Vector ID	VB180313-1293xca						
Vector Name	pPB[Exp]-Puro-EFS>hEEF2K[NM_013302.3]						
Date Created (Pacific Time)	2018-03-13						
Size	8289 bp						
Vector Type	PiggyBac transposon gene expression vector						
Inserted Promoter	EFS						
Inserted ORF	hed ORF hEEF2K[NM_013302.3]						
Inserted Marker	Puro						
Plasmid Copy Number	High						
Antibiotic Resistance	Ampicillin						
Cloning Host	Stbl3 (or alternative strain)						

Vector Map



Vector Components

Name	Position	Size (bp)	Туре	Description	Application notes
5' ITR	■ 1-313	313	ITR piggyBac 5' inverted by		Recognized by PBase transposase; DNA flanked by piggyBac 5' ITR and 3' ITR can be transposed by PBase into TTAA sites.
EFS	■ 337-568	232	Promoter	Human eukaryotic translation elongation factor 1 α1 short form	Strong promoter.
Kozak	■ 593-598	6	Miscellaneous	Kozak translation initiation sequence	Facilitates translation initiation of ATG start codon downstream of the Kozak sequence.
hEEF2K[NM_013302.3]	■ 599-2776	2178	ORF	None	None
rBG pA	■ 2804-3325	522	PolyA_signal	Rabbit beta-globin polyadenylation signal	Allows transcription termination and polyadenylation of mRNA transcribed by Pol II RNA polymerase.

Vector # VB180313-1293xca Page 1 / 4



Name	Position	Size (bp)	Туре	Description	Application notes
CMV promoter	■ 3351-3938	588	Promoter	Human cytomegalovirus immediate early enhancer/promoter	Strong promoter; may have variable strength in some cell types.
Puro	■ 3970-4569	600	ORF	Puromycin resistance gene	Allows cells to be resistant to puromycin.
BGH pA	■ 4613-4837	225	PolyA_signal	Bovine growth hormone polyadenylation signal	Allows transcription termination and polyadenylation of mRNA transcribed by Pol II RNA polymerase.
3' ITR	■ complement (5019-5253)	235	ITR	piggyBac 3' inverted terminal repeat	Recognized by PBase transposase; DNA flanked by piggyBac 5' ITR and 3' ITR can be transposed by PBase into TTAA sites.
Ampicillin	■ 6085-6945	861	ORF	Ampicillin resistance gene	Allows E. coli to be resistant to ampicillin.
pUC ori	■ 7116-7704	589	Rep_origin	pUC origin of replication	Facilitates plasmid replication in E. coli; regulates high-copy plasmid number (500-700).

Note: Components added by user are listed in **bold red** text.

Vector Sequence

1	TTAACCCTAG	AAAGATAGTC	TGCGTAAAAT	TGACGCATGC	ATTCTTGAAA	TATTGCTCTC	TCTTTCTAAA	TAGCGCGAAT	CCGTCGCTGT	GCATTTAGGA
101	CATCTCAGTC	GCCGCTTGGA	GCTCCCGTGA	GGCGTGCTTG	TCAATGCGGT	AAGTGTCACT	GATTTTGAAC	TATAACGACC	GCGTGAGTCA	AAATGACGCA
201	TGATTATCTT	TTACGTGACT	TTTAAGATTT	AACTCATACG	ATAATTATAT	TGTTATTTCA	TGTTCTACTT	ACGTGATAAC	TTATTATATA	TATATTTTCT
301	TGTTATAGAT	ATCATCAACT	TTGTATAGAA	AAGTTGGGCT	CCGGTGCCCG	TCAGTGGGCA	GAGCGCACAT	CGCCCACAGT	CCCCGAGAAG	TTGGGGGGAG
401	GGGTCGGCAA	TTGATCCGGT	GCCTAGAGAA	GGTGGCGCGG	GGTAAACTGG	GAAAGTGATG	TCGTGTACTG	GCTCCGCCTT	TTTCCCGAGG	GTGGGGGAGA
501	ACCGTATATA	AGTGCAGTAG	TCGCCGTGAA	CGTTCTTTTT	CGCAACGGGT	TTGCCGCCAG	AACACAGGCA	AGTTTGTACA	AAAAAGCAGG	CTGCCACCAT
601	GGCAGACGAA	GATCTCATCT	TCCGCCTGGA	AGGCGTTGAT	GGCGGCCAGT	CCCCCGAGC	TGGCCATGAT	GGTGATTCTG	ATGGGGACAG	CGACGATGAG
701	GAAGGTTACT	TCATCTGCCC	CATCACGGAT	GACCCAAGCT	CGAACCAGAA	TGTCAATTCC	AAGGTTAATA	AGTACTACAG	CAACCTAACA	AAAAGTGAGC
801	GGTATAGCTC	CAGCGGGTCC	CCGGCAAACT	CCTTCCACTT	CAAGGAAGCC	TGGAAGCACG	CAATCCAGAA	GGCCAAGCAC	ATGCCCGACC	CCTGGGCTGA
901	GTTCCACCTG	GAAGATATTG	CCACCGAACG	TGCTACTCGA	CACAGGTACA	ACGCCGTCAC	CGGGGAATGG	CTGGATGATG	AAGTTCTGAT	CAAGATGGCA
1001	TCTCAGCCCT	TCGGCCGAGG	AGCAATGAGG	GAGTGCTTCC	GGACGAAGAA	GCTCTCCAAC	TTCTTGCATG	CCCAGCAGTG	GAAGGGCGCC	TCCAACTACG
1101	TGGCGAAGCG	CTACATCGAG	CCCGTAGACC	GGGATGTGTA	CTTTGAGGAC	GTGCGTCTAC	AGATGGAGGC	CAAGCTCTGG	GGGGAGGAGT	ATAATCGGCA
1201	CAAGCCCCCC	AAGCAGGTGG	ACATCATGCA	GATGTGCATC	ATCGAGCTGA	AGGACAGACC	GGGCAAGCCC	CTCTTCCACC	TGGAGCACTA	CATCGAGGGC
1301	AAGTACATCA	AGTACAACTC	CAACTCTGGC	TTTGTCCGCG	ATGACAACAT	CCGCCTGACG	CCGCAGGCCT	TCAGCCACTT	CACTTTTGAG	CGTTCCGGCC
1401	ATCAGCTGAT	AGTGGTGGAC	ATCCAGGGAG	TTGGGGATCT	CTACACTGAC	CCACAGATCC	ACACGGAGAC	GGGCACTGAC	TTTGGAGACG	GCAACCTAGG
1501	TGTCCGCGGG	ATGGCGCTCT	TCTTCTACTC	TCATGCCTGC	AACCGGATTT	GCGAGAGCAT	GGGCCTTGCT	CCCTTTGACC	TCTCGCCCCG	GGAGAGGGAT
1601	GCAGTGAATC	AGAACACCAA	GCTGCTGCAA	TCAGCCAAGA	CCATCTTGAG	AGGAACAGAG	GAAAAATGTG	GGAGCCCCCG	AGTAAGGACC	CTCTCTGGGA
1701	GCCGGCCACC	CCTGCTCCGT	CCCCTTTCAG	AGAACTCTGG	AGACGAGAAC	ATGAGCGACG	TGACCTTCGA	CTCTCTCCCT	TCTTCCCCAT	CTTCGGCCAC
1801	ACCACACAGC	CAGAAGCTAG	ACCACCTCCA	TTGGCCAGTG	TTCAGTGACC	TCGATAACAT	GGCATCCAGA	GACCATGATC	ATCTAGACAA	CCACCGGGAG
1901	TCTGAGAATA	GTGGGGACAG	CGGATACCCC	AGTGAGAAGC	GGGGTGAGCT	GGATGACCCT	GAGCCCCGAG	AACATGGCCA	CTCATACAGT	AATCGGAAGT
2001	ACGAGTCTGA	CGAAGACAGC	CTGGGCAGCT	CTGGACGGGT	ATGTGTAGAG	AAGTGGAATC	TCCTCAACTC	CTCCCGCCTC	CACCTGCCGA	GGGCTTCGGC
2101	CGTGGCCCTG	GAAGTGCAAA	GGCTTAATGC	TCTGGACCTC	GAAAAGAAAA	TCGGGAAGTC	CATTTTGGGG	AAGGTCCATC	TGGCCATGGT	GCGCTACCAC
2201	GAGGGTGGGC	GCTTCTGCGA	GAAGGGCGAG	GAGTGGGACC	AGGAGTCGGC	TGTCTTCCAC	CTGGAGCACG	CAGCCAACCT	GGGCGAGCTG	GAGGCCATCG
2301	TGGGCCTGGG	ACTCATGTAC	TCGCAGTTGC	CTCATCACAT	CCTAGCCGAT	GTCTCTCTGA	AGGAGACAGA	AGAGAACAAA	ACCAAAGGAT	TTGATTACTT
2401	ACTAAAGGCC	GCTGAAGCTG	GCGACAGGCA	GTCCATGATC	CTAGTGGCGC	GAGCTTTTGA	CTCTGGCCAG	AACCTCAGCC	CGGACAGGTG	CCAAGACTGG
2501	CTAGAGGCCC	TGCACTGGTA	CAACACTGCC	CTGGAGATGA	CGGACTGTGA	TGAGGGCGGT	GAGTACGACG	GAATGCAGGA	CGAGCCCCGG	TACATGATGC
2601	TGGCCAGGGA	GGCCGAGATG	CTGTTCACAG	GAGGCTACGG	GCTGGAGAAG	GACCCGCAGA	GATCAGGGGA	CTTGTATACC	CAGGCAGCAG	AGGCAGCGAT
2701	GGAAGCCATG	AAGGGCCGAC	TGGCCAACCA	GTACTACCAA	AAGGCTGAAG	AGGCCTGGGC	CCAGATGGAG	GAGTAAACCC	AGCTTTCTTG	TACAAAGTGG
2801	TGATCCTCAG	GTGCAGGCTG	CCTATCAGAA	$\underline{\mathtt{GGTGGTGGCT}}$	GGTGTGGCCA	ATGCCCTGGC	TCACAAATAC	$\underline{\mathtt{CACTGAGATC}}$	$\underline{\mathtt{TTTTTCCCTC}}$	TGCCAAAAAT
2901	TATGGGGACA	TCATGAAGCC	CCTTGAGCAT	$\underline{\mathtt{CTGACTTCTG}}$	GCTAATAAAG	$\underline{\texttt{GAAATTTATT}}$	TTCATTGCAA	$\underline{\texttt{TAGTGTGTTG}}$	GAATTTTTTG	TGTCTCTCAC
3001	TCGGAAGGAC	ATATGGGAGG	GCAAATCATT	TAAAACATCA	GAATGAGTAT	TTGGTTTAGA	GTTTGGCAAC	ATATGCCCAT	$\underline{\mathtt{ATGCTGGCTG}}$	CCATGAACAA
3101	$\underline{\mathtt{AGGTTGGCTA}}$	TAAAGAGGTC	ATCAGTATAT	GAAACAGCCC	$\underline{\mathtt{CCTGCTGTCC}}$	$\underline{\mathtt{ATTCCTTATT}}$	$\underline{\mathtt{CCATAGAAAA}}$	$\underline{\mathtt{GCCTTGACTT}}$	$\underline{\mathtt{GAGGTTAGAT}}$	TTTTTTTATA
3201	$\underline{\mathtt{TTTTGTTTTG}}$	TGTTATTTTT	$\underline{\mathtt{TTCTTTAACA}}$	$\underline{\mathtt{TCCCTAAAAT}}$	$\underline{\mathtt{TTTCCTTACA}}$	$\underline{\text{TGTTTTACTA}}$	$\underline{\mathtt{GCCAGATTTT}}$	$\underline{\texttt{TCCTCCTCTC}}$	$\underline{\mathtt{CTGACTACTC}}$	CCAGTCATAG
3301	$\underline{\mathtt{CTGTCCCTCT}}$	TCTCTTATGG	AGATCCCTCG	ACCTGCAGCC	CAAGCTTCGC	GTTGACATTG	$\underline{\mathtt{ATTATTGACT}}$	$\underline{AGTTATTAAT}$	$\underline{\mathtt{AGTAATCAAT}}$	TACGGGGTCA

Vector # VB180313-1293xca Page 2 / 4



	mma	000000000000000000000000000000000000000		OFFER CA FRA A C	 magaagaaga		0022000	0000003 EE	20000222022
3401			GGAGTTCCGC		 				
3501			ACGCCAATAG						
3601			CTATTGACGT		 				
3701			CGCTATTACC		 				
3801			GGGAGTTTGT		 				
3901			TCTATATAAG						
4001			AGGGCCGTAC		 				
4101					 TCGGCAAGGT				
4201			GGGGGGGGTG		 				
4301			CCCAAGGAGC		 				
4401			CCGAGCGCGC		 				
4501			GCCCGAAGGA		 				
4601			TCTAGTTGCC		 				
4701			CATCGCATTG		 				
4801			TGCGGTGGGC						
4901			TGTATAGGTT						
5001			TGTTACTTTA		 	-			
5101			AATATAATAA		 	-			
5201			TACGTCACAA		 				
5301			AGTGAGTCGT						
5401			TCGCCAGCTG						
5501			AAGCGCGGCG						
5601			TTCGCCGGCT						
5701			ATGGTTCACG						
5801			AACACTCAAC						
5901			GCGAATTTTA						
6001			TATGTATCCG						
6101			AACTGGATCT		 				
6201			ATCCCGTATT		 				
6301			GATGGCATGA		 				
6401			TAACCGCTTT		 				
6501			GATGCCTGTA		 				
6601			AAAGTTGCAG		 				
6701 6801			CACTGGGGCC		 				
			TGCCTCACTG		 				
6901 7001			AAGATCCTTT						
7101			ATCCTTTTTT						
7201			GAAGGTAACT		 				
7301			TACCTCGCTC		 				
7401			AGCGGTCGGG		 				
7501			CACGCTTCCC		 				
7601			CTTTATAGTC		 				
7701			CCTTTTTACG		 				
7801			AGCTGATACC						
7901			TGGCCGATTC						
8001	TTAGCTCACT		CCCAGGCTTT						
8101			GCTCGAAATT						
8201			TGTGGTTTTT						
0 - 0 1									

Validation by Restriction Enzyme Digestion

Vector # VB180313-1293xca Page 3 / 4



Cutters	Locations	Fragments (bp)		
NaeI	3923, 4366			
NdeI	3012, 3072, 3080, 3604	60, 8, 524, 7697		
ApaLI	6201, 7447	1246, 7043		
ApaLI+NaeI	1704, 5627, 6201, 7447	3923, 574, 1246, 2546		
ApaLI+NdeI	3012, 3072, 3080, 3604, 6201, 7447	60, 8, 524, 2597, 1246, 3854		

Vector # VB180313-1293xca Page 4 / 4